

Portable, Low-Power, Mechanically-Cooled Ge Spectrometer

J.A. Becker¹, C.P. Cork², L. Fabris², and N.W. Madden²

¹*Lawrence Livermore National Laboratory, Livermore, California 94551*

²*Lawrence Berkeley National Laboratory, Berkeley, California 94720*

Field use of Germanium gamma-ray spectrometers is often limited by the availability of liquid nitrogen required for cooling, or by the weight and power of most mechanical coolers. We have constructed a small, portable Germanium spectrometer which takes advantage of a commercial mechanical microcooler[1]. A 5 cm × 5 cm coaxial Germanium detector is cooled to the operating temperature required to achieve the excellent resolution and line shape characteristic of Germanium spectrometers. The Germanium detector is hermetically encapsulated[2] within a thin aluminum can with 2 + atmospheres of nitrogen blanket gas. Two other enabling technologies include the use of an external metal vacuum getter, and a miniature ion-pump both of which are used to attain and maintain the utility vacuum. The result is a Germanium gamma ray spectrometer that has near laboratory resolution, is light in weight and therefore portable, and has low power requirements. Companion low-power electronics were also developed for pulse shaping, histogram collection, and control. Features of the spectrometer are:

- Ge detector
 - 5 cm diameter x 5 cm long
 - Resolution ~ 3.5 keV at $E_\gamma = 662$ keV
 - ~25% efficient relative to 3 x 3" NaI
 - Hermetically encapsulated, nitrogen blanket
- Mechanical cooler, Stirling Cycle
 - Hymatic Engineering
 - Lifetime > 30,000 hours
- Continuous unattended operation
 - ≥ 6 months
- BGO suppression shield (Removable)
 - Improved sensitivity and directionality
 - Reduced background
- Low weight – total = ~ 20 lbs
 - Ge detector + cooler, ~ 10 lbs
 - Detachable suppression shield, 10 lbs
- Low system power
 - 15 Watts, or less
- Integrated electronics
 - 12 bit Nuclear quality ADC
 - 16 Channels data storage (i.e., 16 histograms)
- Operating temperature:
 - - 40 ° C – + 40° C

The concept and present hardware are illustrated in Figure 1. Clearly, reasonably sized Germanium spectrometers are now available for field use. Future developments include installation of a cooled FET preamplifier, in which case we expect laboratory resolution. We anticipate that larger detectors, perhaps 50—80% efficient, can also be cooled with this technique.




Figure 1. Concept of the encapsulated Germanium detector Cryo 3/10 (a), and actualization Cryo 3/25 (b).

References

- 1) The Hymatic Engineering Company Limited.
- 2) N. W. Madden, D. F. Malone, R. H. Pehl, C. P. Cork, P. N. Luke, D. A. Landis, M. J. Pollard, "Germanium Detector Encapsulation", LBL-29857, and Proceedings of International Symposium on Gamma-Ray Line Astrophysics, Paris-Saclay, December 10-13, 1990.

This work was funded in part by DOE/NA, and performed under the auspices of the U.S. Department of Energy by the University of California, Lawrence Livermore National Laboratory and Lawrence Berkeley National Laboratory, Contract No. W-7405-Eng-48 (LLNL), and Contract No. AC03 76SF00098(LBNL).